# JAST (Journal of Animal Science and Technology) TITLE PAGE

# 

ARTICLE INFORMATION	Fill in information in each box below
Article Type	Review article
Article Title	Principal protocols for the processing of cultured meat
Running Title	Processing of cultured meat
Author	Seung Yun Lee <sup>1</sup> , Hea Jin Kang <sup>1</sup> , Da Young Lee <sup>1</sup> , Ji Hyeop Kang <sup>1</sup> , Sivasubramanian Ramani <sup>2</sup> , Sungkwon Park <sup>2</sup> and Sun Jin Hur <sup>1, *</sup>
Affiliation	<sup>1</sup> Department of Animal Science and Technology, Chung- Ang University, Anseong 17546, Korea
	<sup>2</sup> Department of Food Science & Biotechnology, Sejong University, Seoul 05006, Korea
ORCID	Seung Yun Lee
	https://orcid.org/0000-0002-8861-6517
	Hea Jin Kang
	https://orcid.org/0000-0001-6765-3434
	Da Young Lee
	https://orcid.org/0000-0002-3172-0815
	Ji Hyeop Kang
	https://orcid.org/0000-0002-8389-9597
	Sivasubramanian Ramani
	https://orcid.org/0000-0001-9370-6552
	Sungkwon Park
	https://orcid.org/0000-0002-7684-9719
	Sun Jin Hur
	https://orcid.org/0000-0001-9386-5852
Conflicts of interest	The authors declare no potential conflict of interest.

## Upload this completed form to website with submission

Acknowledgements	This research was supported by the Chung-Ang University.
Author contributions	Investigation: Lee SY, Kang HJ, Lee DY, Kang JH, Hur SJ Writing - original draft: Lee SY, Hur SJ Writing - review & editing: Lee SY, Ramani S, Park SK, Hur SJ.
Ethics approval (IRB/IACUC)	This manuscript does not require IRB/IACUC approval because there are no human and animal participants.

## 5 CORRESPONDING AUTHOR CONTACT INFORMATION

For the <u>corresponding</u> author (responsible for correspondence, proofreading, and reprints)	Fill in information in each box below
First name, middle initial, last name	Sun Jin Hur
Email address – this is where your proofs will be sent	hursj@cau.ac.kr
Secondary Email address	
Postal address	Department of Animal Science and Technology, Chung- Ang University, Anseong 17546, Korea.
Cell phone number	+82 10 5441 9438
Office phone number	Tel.: + 82 31 670 4673
Fax number	Fax: + 82 31 670 3108

### 7 Abstract

8 The purpose of this study was to establish a basic principal procedure for the processing of cultured meat. The first 9 stage involved isolating satellite cells from the desired muscle of an animal using enzymatic digestion (i.e., by using 10 proteases, collagenases, and pronases). The second stage involved culturing the isolated muscle satellite cells in a growth medium containing fetal bovine serum and penicillin/streptomycin with growth factors for an optimal period 11 12 of time. The second stage involved a basic method for the isolated muscle cells to proliferate while sub-culturing to 13 further induce differentiation in gelatin-coated culture dishes with the general culture medium. The third stage 14 involved the induction of differentiation of muscle satellite cells or formation of myotubes using myogenic medium. 15 Lastly, the fourth stage involved the identification of cell differentiation or myotube formation (myogenesis) using 16 fluorescent dyes. Moreover, the principle of these protocols can be applied to perform primary culture of animal cells. 17 This study will assist beginners with the technical aspects of culturing meat (isolation, cultivation, and differentiation 18 of muscle satellite cells as well as identification of myotube formation for myogenesis).

- 19
- 20 Key words: Cultured meat, Muscle satellite cell, Cell proliferation, Differentiation, Cell medium



#### 21 INTRODUCTION

22 The negative perception of the livestock industry is increasing due to animal welfare and environmental concerns due 23 to the increase in greenhouse gases and human population. Various alternatives to meat as sources of proteins, such 24 as artificial meat derived from vegetables and the consumption of insects, have been suggested to solve these problems. 25 However, the production of plant-based artificial meat is highly limited in use due to differences in its taste and texture 26 compared to meat. Recently, the production of meat using cultured meat technology has emerged as an alternative to 27 this problem. Cultured meat refers to edible meat obtained through cell proliferation without going through the 28 process of raising livestock and is called in vitro meat or lab-grown meat. This concept of cultured meat has long 29 been suggested, and stem cell technology has been actively applied in this area since its development. Founded in 30 2013 by Professor Mark Post of Maastricht University in the Netherlands, Mosa Meat succeeded in producing beef 31 by cultivating stem cells; thus, the commercial use of cultured meat began in earnest. The use of cultured meat 32 technology not only reduces the energy consumption for meat production to about 45% of that of existing methods 33 of livestock breeding for meat but also reduces greenhouse gas emissions by more than 78-96%, use of land by more 34 than 99%, and use of water by 82-96% compared to the use of animal livestock [1]. Moreover, it has the advantages 35 of preventing or managing animal welfare issues and animal-derived diseases. However, there remain many obstacles 36 to overcome to produce cultured meat for widespread consumption. In particular, since most of the growth factors and culture media used to produce cultured meat are derived from animals, cultured meat does not completely replace 37 38 the use of livestock, and little is known regarding the safety of cultured meat consumption. Besides, it is difficult to 39 clearly resolve this argue that the cultured meat might save the unethically slaughtered livestock since the producing 40 cultured meat need animal sacrifice, suffering and death [2]. Another argue associated with producing cultured meat 41 is that it is unethical to inform unhealthy foods even if it believe that the cultured meat will ethically be produce in 42 the future [3]. In addition, the production of cultured meat is markedly more expensive than the use of conventional 43 livestock; thus, there is a clear limit to its commercial use. Therefore, various attempts are currently being made to 44 develop technologies that produce safe cultured meat inexpensively. At present, various laboratories worldwide are 45 conducting research into cultured meat production, but only a few basic experimental methods have been reported 46 regarding the production of cultured meat. Therefore, this study presents a basic experimental method for the 47 production of cultured meat for commercial use and research (Figure 1). This study is expected to assist with the 48 production of cheap and safe cultured meat.

#### 49 MATERIALS AND METHODS

#### 50 Materials

Analytical grade isoflurane, ethanol (EtOH), bovine serum albumin (BSA), Triton X-100, paraformaldehyde, sucrose, sodium chloride, potassium chloride, sodium phosphate dibasic, potassium phosphate monobasic, 2-methylbutane, protease, and collagenase were purchased from Sigma-Aldrich (St Louis, MO, USA). Dulbeccos' Modified Eagel's Medium (DMEM), fetal bovine serum (FBS), horse serum (HS), and penicillin-streptomycin were purchased from Gibco (Thermo Fisher Scientific, Waltham, MA, USA). Anti-rabbit IgG (H+L) Alexa Fluor 488 conjugate was purchased from Santa Cruz Biotechnology (Santa Cruz, CA, USA), and 4',6-diamidino-2-phenylindole and dihydrochloride were purchased from Invitrogen (Thermo Fisher Scientific, Waltham, MA, USA).

58

### 59 **Preparation of reagents**

60 The cultivating medium for muscle cells was prepared with Ham's F-10, fetal bovine serum (FBS), and antibiotics 61 (penicillin and streptomycin; P/S). Briefly, the final medium was prepared with 100 unit/mL of 20% FBS and 100 62 µg/mL of 1% 1X P/S and stored at 4°C. Myogenic medium for inducing differentiation was prepared with DMEM-63 HG, HS, and antibiotics. Briefly, the final medium was prepared using 2% HS and 1% 1X P/S in DMEM-HG, and 64 stored at 4°C until the end of use. The enzymes used for tissue separation were prepared by adding 3 units of 65 collagenase D along with 4.8 units of dispase II. Finally, to achieve maximum collagenase activity, the digesting 66 enzyme solutions were prepared with DPBS and filtered using a 0.22 µm filter (Whatman<sup>™</sup>, 9913-2502). These 67 solutions were freshly prepared each time prior to digestion. Calcium chloride (1.47 g) was added to 100 mL of 68 sterilized distilled water to form a 100 mM concentration, which was then filtered using a 0.22 µm filter and stored. 69 The 10 ng/mL basic fibroblast growth factor (bFGF) solution was prepared with 10 ng of bFGF by adding 1 mL of 70 0.1% BSA to PBS. A total of 100 mL of 75% EtOH was prepared using 100% EtOH (75 mL) and sterilized ultra-71 pure water (25 mL). Sterile single-use disposable plastic wares, such as cell strainers (40, 70, and 100 µm), petri 72 dishes, polypropylene centrifuge tubes (15 and 50 mL), and T-flasks (T75), were used. A clean environment was 73 maintained that included a water bath, centrifuge, incubator, and microscope.

74

#### 75 Isolation of muscle cell population

- The institute is an imal care guidelines, and clean the surface of desired
   muscle region with EtOH spray. Then, resect the muscle (from the hind limbs) of the animal using surgical tools.
- 78 <u>2. While collecting the muscle, take care to ensure that hair and any other external materials are not included.</u>

- $\frac{3. \text{ Sterilize the collected (approximately 10 g) muscle tissue (approximately <math>1 \times 1 \text{ cm}$ ) with 75% EtOH and wash
- 80 three times with ice cold 1% P/S in DPBS (v/v).
- 4. Remove the fat and other binding tissues, and mince the muscle tissues (to approximately 2-5 g and 5 × 5 mm)
   using sterilized scissors or a razor blade (35 mm).
- 5. To isolate muscle cells, use a 1:3 ratio of the prepared digestive enzyme solution for the digestion of the external
   intercellular matrix by combining the cells and solution in a 15 mL tube at 37°C for 30-60 min.
- 85 <u>\* To prepare digestive enzymes, mix the digestive enzymes following 1 mL of 3 U collagenase D (final concentration: 1.5 U collagenase D), 1 mL of 4.8 U dispase II (final concentration: 2.4 U dispase II) and 50 μL
  </u>
- 87 of 100 mM calcium chloride (final concentration: 2.5 mM calcium chloride).
- 88 Generally, several types of enzymes can be used to digest muscle tissues, such as proteases, collagenases (e.g.,
- 89 <u>collagenase VI</u>), trypsin, and pronase, which can affect the extracellular matrix.
- 90 <u>6. After enzymatic digestion, remove non-singular mass using a strainer of pore sizes 100, 70, and 40 μm.</u>
- 91 <u>7. Recover the filtered cells by differential centrifugation ( $350 \times g$  for 5 min at 4°C). Discard the supernatant and</u>
- 92 repeat this step twice to wash off enzymes and condition the cells by further digestion and pH alterations using
   93 culture medium with 20% FBS.
- 8. Recover the pelleted cells (approximately 1 g) and suspend as 10 mL of culture medium (20% FBS in Ham's F10 with 1% P/S) and 10 ng/mL bFGF in a T75 flask.
- 96 <u>\* bFGF promotes satellite cell proliferation and inhibits the differentiation of myoblasts [4]. In addition, some</u>
- 97 studies have suggested that the addition of leukemia inhibitory factor, transforming growth factor- $\beta$ , and insulin-
- 98 like growth factor-I can improve animal muscle satellite [5, 6].
- 99 <u>9. Incubate the suspended cells in the T75 flask at 37°C in a 5% CO<sub>2</sub> incubator for 2 h.</u>
- 100 Therefore, this step can be applied to isolate whole muscle cell populations with many satellite cells from animals101 (Figure 2).
- 102 Further steps were carried out to culture the isolated muscle cells to produce cultured meat.
- 103

#### 104 Culture of muscle satellite cells

- 105 <u>1. Transfer the supernatant obtained following the above steps (isolation of muscle cells) in the incubated T-75</u>
   106 <u>flask into either a 0.1% gelatin- or Matrigel<sup>©</sup>-coated T75 flask.</u>
- 107 <u>2. Collect the remaining cultured cells (No. 9) in the T75 flask using culture medium (approximately 3-5 mL) twice</u>
- 108and then transfer into a new T75 flask.

109 <u>3. Incubate the transferred cells in the gelatin- or Matrigel-coated T75 flasks without medium change at 37°C in a</u>

110  $\underline{5\% \text{ CO}_2 \text{ incubator for 72 h.}}$ 

- 111 <u>4. After 72 h, wash the incubated cells and change the culture medium at intervals of two to four days.</u>
- 5. For subculturing, remove the old medium and wash with DPBS twice and add 1-2 mL of 0.05% trypsin, followed
   by incubation at 37°C in a 5% CO<sub>2</sub> incubator for 2-3 min.
- 114 <u>6. After incubation, to maintain cell viability, add the 0.5 mL of DPBS to neutralize trypsin and stop the digestion</u>
- 115 process. From the neutralized trypsin solution, separate the cells using a centrifuge at  $300 \times g$  for 5 min at 4°C.
- 116 <u>7. After centrifugation, suspend the pellet in 7-10 mL of culture medium (20% FBS in Ham's F-10 with 1% P/S).</u>
- 8. Incubate the suspended cells in 0.1% gelatin-coated T75 flasks for 1 h, repeating twice to separate only the fiber
   cells and sample muscle satellite cells.
- \* While fibroblasts and skeletal muscle satellite cells were cultured in 0.1% gelatin-coated T75 flasks due to
   different culture characteristics, fibroblasts quickly attached to the bottom of the cell culture flask and muscle
- satellite cells were slow to attach and remained in the suspended fluid [7].
- 9. Based on the characteristics of muscle satellite cells, collect the supernatant as satellite cells (approximately 7 10 mL) in the T75 flask and transfer into a 0.1% gelatin-coated T75 flask.
- 124 <u>10. This step was repeated two to three times using culture medium (20% FBS in Ham's F-10 with 1% P/S) to</u>
   125 <u>collect the remaining cells, which might have been rich in satellite cells.</u>
- 126 <u>11. Seed the isolated satellite cells at a density of 0.3 × 10<sup>6</sup> cells with 7-10 mL of culture medium (20% FBS in
   127 Ham's F-10 with 1% P/S and 10 ng/mL bFGF).
  </u>
- 128 <u>12. Incubate the cells to grow at  $37^{\circ}$ C in a 5% CO<sub>2</sub> incubator.</u>
- 129 <u>13. Maintain the cells by subculturing until that passage number of cells is 13-15 (Figure 3).</u>

- 131 Inducing muscle satellite cell differentiation
- 132 <u>1. To induce the differentiation of muscle satellite cells, change the culture medium to myogenic medium (2% HS</u>
- 133 <u>and 1% P/S).</u>
- 134 <u>2. Aspirate carefully the cultivation medium, and add the cultured muscle satellite cells to 10 mL of myogenic</u>
   135 medium (2% HS and 1% P/S) and then incubate at 37°C and 5% CO<sub>2</sub> for one week.
- 3. After one week, observe the differentiation of cells and formation of myotubes using fluorescent dye reagents
   (Figure 4).
- 138

- 139 Identification of differentiated satellite cells/myogenesis using biomarkers
- The biomarkers of satellite and differentiated cells were used to identify the satellite and differentiated cells as well as myotube formation as indicators of myogenesis. Here, the differentiation cells <u>were</u> quantified using the cell fusion index, which refers to the percentage of nuclei inside myotubes compared to the total count of nuclei in the sample
- 143 [6]. Therefore, the fusion index was assessed by determining the visualization of nuclear staining, such as by using
- 144 fluorescent dyes (e.g., bisBenzimide H 33342 trihydrochloride and Mito Tracker Red CMX-Ros) [8] and Giemsa
- 145 staining.
- 146 <u>1. Prepare the fluorescent dye solution (200 μM of Mito Tracker Red CMX-Ros) from the stock right before</u>
   staining the cells and immediately use as a fresh preparation to avoid the loss of fluorescence ability.
- 148 <u>\* prepare a 200 μM fluorescent dye stock solution by dissolving 50 μg of Mito Tracker Red CMX-Ros dye in 470</u>


149  $\mu$ L of DMSO, which is then stored at -20°C.

- 150 <u>2. Treat the differentiated cells in the T75 flasks with 10 mL of medium with labeling solution (2.5 μL of 200 μM</u>
- Mito Tracker Red CMX-Ros (final concentration: 50 nM) and 10 μL of Hoechst 33342 dye (final concentration:
   1 μM)) for 30 min at 37°C.
- <u>3. After incubation, remove the labeling solution and Hoechst 33342 dye, and wash the cells using cell medium</u>
   <u>twice.</u>
- 155 <u>4. Fix the Labeled cells with 4% formaldehyde for 15 min at 37°C, followed by washes in DPBS.</u>
- 156 <u>5. Observe differentiated cells or myotubes using fluorescence microscopy (Figure 4).</u>
- 157
- 158 References
- 159 1. Tuomisto, H. L., & Teixeira de Mattos, M. J. (2011). Environmental impacts of cultured meat production.

160 Environmental Science & Technology, 45(14), 6117-6123.


- 161 2. Schaefer, G. O., & Savulescu, J. (2014). The Ethics of Producing In Vitro Meat. Journal of applied philosophy,
- 162 31(2), 188–202.
- 163 3. Bhat, Z. F., Morton, J. D., Mason, S. L., Bekhit, A. E. D. A., & Bhat, H. F. (2019). Technological, regulatory,
- and ethical aspects of *in vitro* meat: A future slaughter-free harvest. Comprehensive Reviews in Food Science
- 165 and Food Safety, 18(4), 1192-1208.
- 166 4. Guthridge M, Wilson M, Cowling J, Bertolini J, Hearn MT. The role of basic fibroblast growth factor in skeletal
- 167 muscle regeneration. Growth Factors 1992. 6.

- 168 5. Allen RE, Boxhorn LK. Regulation of skeletal muscle satellite cell proliferation and differentiation by
- transforming growth factor-beta, insulin-like growth factor I, and fibroblast growth factor. J. Cell. Physiol. 1989.
  170 138.
- 171 6. Ben-Arye T, Levenberg S. Tissue engineering for clean meat production. Front. sustain. food syst. 2019. 3.
- 172 7. Hindi L, McMillan JD, Afroze D, Hindi SM, Kumar A. Isolation, culturing, and differentiation of primary
- 173 myoblasts from skeletal muscle of adult mice. Bio Protoc. 2017. 7.
- 174 8. Miyake T, McDermott JC, Gramolini AO. A method for the direct identification of differentiating muscle cells
- by a fluorescent mitochondrial dye. PLoS One. 2011. 6.





178 Fig. 1. Overview of muscle satellite cell culture.





181 Fig. 2. Isolation of muscle satellite cells from animals step.





185 Fig. 4. Differentiation and identification of muscle satellite cells step.

